

# Product Categories

## 1 Fingertip Pulse Oximeter

- 01 AS-301 AS-302
- 02 AS-303 AS-304
- 03 AS-301-L AS-302-L
- 04 AS-304-L

## 2 Fetal Doppler

- 05 AF-703 AF-706
- 06 AF-703-L AF-706-L
- 07 AF-708

## 3 Infrared Thermometer

- 08 AT-601 AT-606
- 09 AT-611


## 4 Compressor Nebulizer

- 10 AN-821 AN-831
- 11 AN-832

## 5 Blood Pressure Monitor

- 12 AB-501 AB-501-B
- 13 AB-503
- 14 AB-511M AB-512L
- 15 AB-513L

## AS-301




- OLED colour screen, bar graph and waveform
- Automatic 4 directions screen rotation
- Beep sound and setting function
- Easy operation & portable design
- Automatically power off
- Accuracy  
SpO<sub>2</sub>(70%~99%): $\pm 2\%$   
PR(30BPM~250BPM): $\pm 1\text{BPM}$  or  $\pm 1\%$
- Low perfusion $\leq 0.2\%$
- Colours: Blue, Light blue, Pink, Purple



## AS-302

- OLED colour screen, bar graph and waveform
- Automatic 4 directions screen rotation
- Beep sound and setting function
- Easy operation & portable design
- Automatically power off
- Accuracy  
SpO<sub>2</sub>(70%~99%): $\pm 2\%$   
PR(30BPM~250BPM): $\pm 1\text{BPM}$  or  $\pm 1\%$
- Low perfusion $\leq 0.2\%$
- Colours: Blue, Light blue, Pink, Purple, Black





G-Sensor



Strong anti-interference



## AS-303

- OLED colour screen, bar graph and waveform
- Automatic 4 directions screen rotation
- Beep sound and setting function
- Easy operation & portable design
- Automatically power off
- Accuracy
  - SpO2(70%~99%):  $\pm 2\%$
  - PR(30BPM~250BPM):  $\pm 1\text{BPM}$  or  $\pm 1\%$
- Low perfusion  $\leq 0.2\%$
- Colours: White

## AS-304



- OLED colour screen, bar graph and waveform
- Automatic 4 directions screen rotation
- Beep sound and setting function
- Easy operation & portable design
- Automatically power off
- Accuracy
  - SpO2(70%~99%):  $\pm 2\%$
  - PR(30BPM~250BPM):  $\pm 1\text{BPM}$  or  $\pm 1\%$
- Low perfusion  $\leq 0.2\%$
- Colours: Purple, Yellow

## AS-301-L



- Red LED screen, bar graph
- Easy operation & portable design
- Automatically power off
- Beep sound
- Accuracy  
SpO<sub>2</sub>(70%~99%): ±2%  
PR(30BPM~250BPM): ±1BPM or ±1%  
Low perfusion ≤0.2%
- Colours: Blue, Light blue, Pink, Purple



## AS-302-L

- Red LED screen, bar graph
- Easy operation & portable design
- Automatically power off
- Beep sound
- Accuracy  
SpO<sub>2</sub>(70%~99%): ±2%  
PR(30BPM~250BPM): ±1BPM or ±1%  
Low perfusion ≤0.2%
- Colours: Blue, Light blue, Pink, Purple, Black








## AS-304-L



- Red LED screen, bar graph
- Easy operation & portable design
- Automatically power off
- Beep sound
- Accuracy
  - SpO<sub>2</sub>(70%~99%): ±2%
  - PR(30BPM~250BPM): ±1BPM or ±1%
  - Low perfusion ≤0.2%
- Colours: Purple, Yellow



## AF-703



- TFT display
- Can detect FHR for pregnancy of above 12 weeks
- Built-in loudspeaker , clear and loud voice
- Rechargeable battery,battery status indicator
- Automatic power off in 1 minute without operation
- Water proof probe
- Headphone socket
- Ultrasound Frequency:2.5MHz
- Battery options:4xAA / chargeable batteries



## AF-706

- TFT display
- Can detect FHR for pregnancy of above 12 weeks
- Built-in loudspeaker , clear and loud voice
- Rechargeable battery,battery status indicator
- Automatic power off in 1 minute without operation
- Water proof probe
- Headphone socket
- Ultrasound Frequency:2.5MHz
- Battery options:4xAA / chargeable batteries





## AF-703-L

- LCD display
- Can detect FHR for pregnancy of above 12 weeks
- Built-in loudspeaker , clear and loud voice
- Rechargeable battery,battery status indicator
- Automatic power off in 1 minute without operation
- Water proof probe
- Headphone socket
- Ultrasound Frequency:2.5MHz
- Battery options:4xAA / chargeable batteries



## AF-706-L

- LCD display
- Can detect FHR for pregnancy of above 12 weeks
- Built-in loudspeaker , clear and loud voice
- Rechargeable battery,battery status indicator
- Automatic power off in 1 minute without operation
- Water proof probe
- Headphone socket
- Ultrasound Frequency:2.5MHz
- Battery options:4xAA / chargeable batteries



## AF-708



- All in one integrated design, small beautiful appearance
- Light in weight, portable to carry, easy operation
- Mainly used to detect the fetal heartbeat rate (FHR) and the sound of the fetal heartbeat.
- The probe taking the principal of fetal ultrasound to collect the fetal heart beat signals of fetus
- Strong ability of processing fetal heart signals
- Auto Shut-off in 1 minute
- Ultrasound Frequency: 2.5MHz



*New*



## AT-601

- Fever alarming
- 9 sets of memories
- Unique mini design
- LCD display screen
- °C/°F free conversion
- Measurement range: 34°C~42.2°C



## AT-606

- Fever alarming
- 9 sets of memories
- LCD display screen
- °C/°F free conversion
- Measurement range: 34°C~42.2°C



# Non-contact Infrared thermometer



## AT-611



- Voice system
- Setting, Fever alarming
- 32 sets of memories
- Tri-color Backlight LCD screen
- °C/°F Temperature units free conversion
- Three temperature modes for Body, Object and Environment

*New*



## AN-821

- Power: 100V@60Hz, 120V@60Hz, 230V@50Hz, 220V@60Hz (Option)
- Extreme Pressure:  $\geq 150$  kpa
- Noise level:  $\leq 65$  dB
- Nebulizing rate:  $\geq 0.2$  mL/min
- Nebulizing Pressure: 86KPa~106KPa
- Air flow:  $\geq 10$  L/min
- MMAD: 2.5 $\mu$ m
- Size: 269mm\*169.5mm\*96mm



## AN-831

- Power: 100V@60Hz, 120V@60Hz, 230V@50Hz, 220V@60Hz (Option)
- Extreme Pressure:  $\geq 150$  kpa
- Noise level:  $\leq 65$  dB
- Nebulizing rate:  $\geq 0.2$  mL/min
- Nebulizing Pressure: 86KPa~106KPa
- Air flow:  $\geq 10$  L/min
- MMAD: 2.5 $\mu$ m
- Size: 164mm\*160mm\*98mm



# Compressor Nebulizer



## AN-832



- Power: 100V@60Hz, 120V@60Hz, 230V @50Hz, 220V@60Hz(Option)
- Extreme Pressure:  $\geq 150\text{kpa}$
- Noise level:  $\leq 65\text{dB}$
- Nebulizing rate:  $\geq 0.2\text{mL/min}$
- Nebulizing Pressure: 86KPa~106KPa
- Air flow:  $\geq 10\text{L/min}$
- MMAD: 2.5 $\mu\text{m}$
- Size: 199mm\*150mm\*95mm

*New*





## AB-501



- Measurement Method: Oscillometric Method
- Strong Anti-Interference Ability
- Display: 3.5inch Color LCD Screen, Touch Screen
- Cuff Size: 22-32cm
- Accuracy
  - Pressure:  $\pm 3$  mmHg
  - Pulse:  $\pm 5\%$
- Memories: 2\*120
- Trend Chart Analysis
- Power Supply: 4× AA Batteries



Touch Screen

## AB-501-B

- Measurement Method: Oscillometric Method
- Strong Anti-Interference Ability
- Display: 3.2 inch Color LCD Screen
- Cuff Size: 22-32cm
- Accuracy
  - Pressure:  $\pm 3$  mmHg
  - Pulse:  $\pm 5\%$
- Memories: 2\*80
- Power Supply: 4× AA Batteries



Back light display



# Blood Pressure Monitor



## AB-503



- Measurement Method : Oscillometric Method
- Strong Anti-Interference Ability
- Display: 2.8 inch LCD Screen
- Cuff Size: 22-32cm
- Accuracy
  - Pressure:  $\pm 3$  mmHg
  - Pulse:  $\pm 5\%$
  - Memories: :2\*80
- Power Supply: 4\* AA Batteries



## AB-511M



- Measurement Method: Oscillometric Method
- Strong Anti-Interference Ability
- Display: 3.9inch white color LCD with big front display
- Cuff Size: 22-32cm
- Accuracy  
Pressure:  $\pm 3$  mmHg  
Pulse:  $\pm 5\%$   
Memories: 2\*60
- Power Supply: 4× AA Batteries
- Display options: 3.2inch(AB-511S)

*New*

## AB-512L



- Measurement Method: Oscillometric Method
- Strong Anti-Interference Ability
- Display: 4.5inch white color LCD with big front display
- Cuff Size: 22-32cm
- Accuracy  
Pressure:  $\pm 3$  mmHg  
Pulse:  $\pm 5\%$   
Memories: 2\*60
- Power Supply: 4× AA Batteries
- Display options: 3.9inch(AB-513M)/3.2inch(AB-513S)

# Blood Pressure Monitor



## AB-513L



- Measurement Method: Oscillometric Method
- Strong Anti-Interference Ability
- Display: 4.5inch white color LCD with big front display
- Cuff Size: 22-32cm
- Accuracy
  - Pressure:  $\pm 3$  mmHg
  - Pulse:  $\pm 5\%$
- Memories: 2\*60
- Power Supply: 4\* AA Batteries
- Display options: 3.9inch(AB-513M)/3.2inch(AB-513S)

*New*



**HAORUN MEDICAL PRODUCTS CO.,LTD**

TEL 0086-574-87706193

WHATSAPP/WECHAT/TEL:008613566542343

EMAIL:SALES@HAORUNMED.COM

WWW.HAORUNMED.COM